

New!

Introducing Gus Gear's LOCK 2000 Vest



Our updated Central Line Vest increases comfort and ease of use, without compromising our unrivaled securement. Here's what to expect:

Central Line Vest



Shoulder Strap Width
0.75 in / 1.9 cm

Front Flap Height
4.25 in / 10.8 cm



Armpit to Vest Bottom
2.75 in / 7 cm



Back Panel Height
6 in / 15.2 cm



Velcro closures in center and flap

LOCK 2000 Vest



Shoulder Strap Width
0.5 in / 1.3 cm

Front Flap Height
3.5 in / 8.9 cm



Armpit to Vest Bottom
1.75 in / 4.4 cm



Back Panel Height
2.75 in / 7 cm



Snap closures in center and flap

Size 3 Vest Comparison

What's New & Why It Matters

Reduced vertical footprint	Enhances breathability and fit
Slimmer shoulder and chest straps	Optimizes patient comfort
Upgraded fabric	Boosts comfort & compliance
Simplified center closure	Ensures consistent results
Unchanged Dual Safe Strip	Provides unrivaled securement
Expanded size range	Fits your entire patient population

Full Size Range Comparison Chart

Key:	Central Line Vest	New LOCK 2000 Vest
------	-------------------	--------------------

	Chest Measurement		Front Flap Height		Shoulder Strap Width		Armpit to Vest Bottom		Back Panel Height	
	in	cm	in	cm	in	cm	in	cm	in	cm
Size 1	13.5 - 16	34 - 40.5	3	7.6	0.75	1.9	1.625	4.1	5.25	13.3
	13 - 16.5	33 - 42	3.125	7.9	0.5	1.3	1.5	3.8	2.75	7
Size 2	16 - 19	40.5 - 48	3.5	8.9	0.75	1.9	1.875	4.8	5.5	14
	16 - 19.5	40.5 - 49.5	3.5	8.9	0.5	1.3	1.75	4.4	2.75	7
Size 3	19 - 22.5	48 - 58	4.25	10.8	0.75	1.9	2.75	7	6	15.2
	18 - 24	45.5 - 61	3.5	8.9	.05	1.3	1.75	4.4	2.75	7
Size 4	23 - 28	58 - 71	5	12.7	1	2.5	3.625	9.2	6.5	16.5
	21.5 - 30	54.5 - 76	3.5	8.9	0.5	1.3	1.75	4.4	2.75	7
Size 5	28 - 34	71 - 86.5	5.75	14.6	1	2.5	3.75	9.5	7	17.8
	28 - 41	71 - 104	4.75	12.1	1	2.5	2.75	7	3.75	9.5
Size 6	34 - 42	86.5 - 107	6.5	16.5	1	2.5	3.75	9.5	7.5	19.1
	37 - 56	94 - 142	4.75	12.1	1	2.5	2.75	7	3.75	9.5